

Assist-a-Lift Capri



POWERLIFT CHAIR

Made in USA

Assists in sitting and standing

3 Way mechanism

Easy to use hand controller

High quality Birch construction

Provides comfortable seat with many rise and recline positions

to suit the user

Available in a range of colours

Dimensions:

Seat width: 510mm

Seat depth: 480mm

Seat height: 460mm

Back height; 690mm

Overall width: 790mm

Maximum user weight: 140kg

\$1,595

HIRE: \$294
FOR THE FIRST 4 WEEKS

Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday