

Assist-a-Lift Supreme

POWERLIFT CHAIR



Made in USA

Powered armchair assists the user in sitting or standing

Provides a comfortable seat with infinite rise and recline positions

Supplied with head covers and magazine pocket

Overfilled seamed back provides increased lumbar support

Hand held controller is simple to use

Rise or recline motion stops at any point when button is released

Chair frame manufactured from American Birch

Lift and recline mechanism

manufactured from high grade steel

Available in a range of colours

Dimensions:

Seat width: 535mm

Seat depth: 545mm

Seat height: 470mm

Back height: 710

Overall width: 815

Maximum user weight: 160kg

\$1,895

\$294* to hire

FOR THE FIRST 4 WEEKS

Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday