

# BATHROOM AIDS

## Shower Chair



Tubular Steel fabrication  
Moulded plastic seat and backrest  
Adjustable Seat Height 410 - 540mm  
Seat area 440(w) x 430(d)mm  
Zinc treated frame to resist rusting  
Armrests to assist transfers

**HIRE:**  
from \$64\* for the first 4 weeks

**PURCHASE:**  
**\$125**

## Shower Stool



Tubular Steel fabrication  
Moulded plastic seat  
Adjustable Seat Height 480 - 610mm  
Seat area 340(w) x 300(d)mm  
Zinc treated frame to resist rusting  
Armrests to assist transfers

**HIRE:**  
from \$64\* for the first 4 weeks

**PURCHASE:**  
**\$125**

## Transfer Bench



Designed for use with two legs inside  
& two legs outside the bath  
Adjustable seat height 440 - 560mm  
Seat area 800(w) x 360(d)mm  
Optional extension legs & suction caps  
Backrest and grab rail for support  
Fits most bath tubs

**HIRE:**  
from \$94\* for the first 4 weeks

**PURCHASE:**  
**\$190**

## Timber Bathboard



Slatted Timber construction  
Chrome grab rail to assist transfers  
Sliding block aids stable fitting  
Fits most bath tubs  
Seat area 710(w) x 310(d)mm  
Optional smaller sizes are available  
to suit thinner bath tubs

**HIRE:**  
from \$64\* for the first 4 weeks

**PURCHASE:**  
**\$125**

\*NB: Terms, conditions and deposits apply to hired goods

## GRAB RAILS:

We stock a range of steel grab rails which may be fitted to stable surfaces in the bathroom (walls & cabinets etc) to give support and aid stability. We supply specialised stainless steel rails, 'made to order' in various shapes to suit your particular requirements.

## ACCESSORIES:

Also available are a wide variety of bathing accessories including;

- Bath & Shower Mats
- Shower Sandle 'Toe & Foot Washers'
- Hand showers (temporary & permanent)
- Tap Turners
- Long handled sponges
- Sitz Baths



**GILLESPIES**  
**HIRE & SALES SERVICE**  
*"Comfort you deserve"*

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

**Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Middy**