

BATHROOM AIDS

Shower Chair



Tubular Steel fabrication
Moulded plastic seat and backrest
Adjustable Seat Height 410 - 540mm
Seat area 440(w) x 430(d)mm
Zinc treated frame to resist rusting
Armrests to assist transfers

HIRE:
from \$50* for the first 4 weeks

PURCHASE:
\$125.00

Shower Stool



Tubular Steel fabrication
Moulded plastic seat
Adjustable Seat Height 480 - 610mm
Seat area 340(w) x 300(d)mm
Zinc treated frame to resist rusting
Armrests to assist transfers

HIRE:
from \$50* for the first 4 weeks

PURCHASE:
\$125.00

Transfer Bench



Designed for use with two legs inside
& two legs outside the bath
Adjustable seat height 440 - 560mm
Seat area 800(w) x 360(d)mm
Optional extension legs & suction caps
Backrest and grab rail for support
Fits most bath tubs

HIRE:
from \$80* for the first 4 weeks

PURCHASE:
\$190.00

Timber Bathboard



Slatted Timber construction
Chrome grab rail to assist transfers
Sliding block aids stable fitting
Fits most bath tubs
Seat area 710(w) x 310(d)mm
Optional smaller sizes are available
to suit thinner bath tubs

HIRE:
from \$50* for the first 4 weeks

PURCHASE:
\$115.00

*NB: hired goods

GRAB RAILS:

We stock a range of steel grab rails which may be fitted to stable surfaces in the bathroom (walls & cabinets etc) to give support and aid stability. We supply specialised stainless steel rails, 'made to order' in various shapes to suit your particular requirements.

ACCESSORIES:

Also available are a wide variety of bathing accessories including;

- Bath & Shower Mats
- Shower Sandle 'Toe & Foot Washers'
- Hand showers (temporary & permanent)
- Tap Turners
- Long handled sponges
- Sitz Baths



PHONE (02) 9411 2180

11-13 Elizabeth Street (PO Box 45) ARTARMON NSW 1570 AUSTRALIA