



**GILLESPIES**  
HIRE & SALES SERVICE

"Comfort you deserve"

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday

*\$490 to hire  
4 weeks*

# CTM HS6200 Folding Power Chair



Code SC0195B

RRP **\$3495**

- Sturdy, Foldable Power Chair
- Connectorless Batteries
- Height Adjustable/ Flip-up Armrests
- Side Exchangeable Controller
- Length Adjustable Controller
- Maximum Speed 8km
- Capacity 127kg

SPECIFICATION	HS 6200
Overall Length	1075mm / 42.3"
Overall Width	650mm / 25.6"
Overall Height	900mm / 35.4"
Wheel Size Front	200mm / 8"
Wheel Size Rear	320mm / 13"
Weight With Batteries	74.4kg / 162lb
Max Speed	8kmph / 5mph
Weight Capacity	127kg / 280lb
Ground Clearance	46 mm / 1.8"
Grade Climbable	10°
Curb Climbing	28mm / 1.1"
Turning Radius	1100mm / 43.3"
Brake	Electro-Mechanical
Seat Width	431.8mm / 17"
Drive Train	2 Motor Rear Wheel Drive
Battery Weight	14.7kg / 32.5lb x 2
Motor Size	4600rpm, 420w
Travel Range	40km / 25m
Battery	12v 36amp x 2
Battery Charger	5amp
Electronics	PG VSI
Seat Type	Padded, Breathable Ballistic Nylon Back Rest

