


Classic FS Bed



'DAYS' Classic FS Bed

FEATURES

- 4-section mattress support, 2-motor adjustment, easy transport
- Client position assistance pole with triangle
- Wooden head end, foot boards, side rails in beech finish
- Handset includes locking facility to maintain position
- Removable mattress retainers at head and foot sections
- Ergonomic adjustment of head and foot sections
- Additional positioning of foot section by ratchet fittings
- Emergency lowering by quick-action bolt mechanism
- Transport lock
- Dual action headrest and knee break adjustment on handset, to help reduce shear when profiling

TECHNICAL DATA

- Total weight....(4-piece) 91.5kg
- Heaviest part....20.1kg
- Outer dimensions W/L....102.3 x 212.5cm
- Max safe working load....175kg
- Max patient weight....140kg
- Height adjustment....40cm to 80cm
- Adjustment angle of back section....0° - 70°
- Adjustment angle of foot section....0° - 20°

QUALITY

- EN 1970:2000 | DIN EN 1970:2000-12
- EN 60601-2-38:1996 + A1:2000 (as basis)
- DIN EN 60601-2-38:2001 (as basis)
- CE 93/42/EEC on medical devices

*Great value,
Fast delivery!*



HIRE: \$540*
FOR THE FIRST 4 WEEKS
Hire cost includes mattress,
side rails and self-help pole

\$2,200

SALE PRICE MATTRESS NOT INCLUDED

*Terms, conditions and deposits apply to hired goods

Code: DAYSFSBEDSULTRA



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday