



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday



The EC-B08 sets a new standard in Luxury Day chair design. The three way segmented Backrest is expertly contoured to provide both Lumbar and Neck support without compromise in overall comfort. The Frame is a sturdy steel construction, zinc treated & powder coated for durability and resistance from corrosion. Legs are adjustable to provide a seat height range from 43cm –53cm. At the push of a button. Seat depth can also be adjusted via screws underneath. Medium density padding of the highest quality has been used in the seat to help eliminate strain on the users lower spinal area.

The backrest reclines to four fixed positions by gently lifting the armrests until the most suitable position is found. Armrests are also upholstered and softly padded for the ultimate in user comfort.

SPECIFICATIONS:

- Easy to clean wipe-over upholstery
- Adjustable height legs (43cm to 53cm)
- Adjustable height backrest (35cm to 48cm)
- Reclining Backrest — 4 positions
- Seat width — 53.5cm
- Seat Depth — 46cm (adjustable)
- SAFE WORKING LOAD 130Kg

\$525

HIRE: \$204 For first 4 weeks

Terms, conditions and deposits apply to hired goods



Recline—4 positions

Gently Lift



DAY CHAIR (EC-B08)

The ultimate in daily living comfort!