


# Classic FS Bed



## 'DAYS' Classic FS Bed

### FEATURES

- 4-section mattress support, 2-motor adjustment, easy transport
- Client position assistance pole with triangle
- Wooden head end, foot boards, side rails in beech finish
- Handset includes locking facility to maintain position
- Removable mattress retainers at head and foot sections
- Ergonomic adjustment of head and foot sections
- Additional positioning of foot section by ratchet fittings
- Emergency lowering by quick-action bolt mechanism
- Transport lock
- Dual action headrest and knee break adjustment on handset, to help reduce shear when profiling

### TECHNICAL DATA

- Total weight....(4-piece) 91.5kg
- Heaviest part....20.1kg
- Outer dimensions W/L....102.3 x 212.5cm
- Max safe working load....175kg
- Max patient weight....140kg
- Height adjustment....40cm to 80cm
- Adjustment angle of back section....0° - 70°
- Adjustment angle of foot section....0° - 20°

### QUALITY

- EN 1970:2000 | DIN EN 1970:2000-12
- EN 60601-2-38:1996 + A1:2000 (as basis)
- DIN EN 60601-2-38:2001 (as basis)
- CE 93/42/EEC on medical devices

*Great value,  
Fast delivery!*



**HIRE: \$540\***

**FOR THE FIRST 4 WEEKS**  
Hire cost includes mattress,  
side rails and self-help pole

**\$2,200**

**SALE PRICE MATTRESS NOT INCLUDED**

\*Terms, conditions and deposits apply to hired goods

Code: DAYSFSBEDSULTRA



**GILLESPIES**  
**HIRE & SALES SERVICE**  
"Comfort you deserve"

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

**Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday**

# HOSPITAL BEDS & PATIENT LIFTERS

We hire and sell a range of bedroom equipment for both temporary & permanent requirements. Remote control; electrically operated hospital beds and patient lifters are available as well as numerous aids to assist transfers between bed and mobility aid, or in the case of a patient lifter, to assist a carer in the safe transfer of a patient with a minimum of strain and discomfort to either party.

## HOSPITAL BEDS

We have a range of hospital beds and mattresses to suit most needs. We also carry a range of pressure care overlays and alternating air mattresses for those requiring additional relief. (Please refer to our 'Pressure Care' brochure) Please feel welcome to call and discuss your specific needs with our experienced staff. Hospital Bed 'safety sides' and 'self help poles' are optionally available.



**HIRE:** \$300\*  
for the first 4 weeks

**PURCHASE:** Days Classic FS \$2,955

## PATIENT LIFTERS

We hire electric patient lifters with slings to aid carers in the safe and convenient transfer of mobility impaired persons. The patient lifters are designed to meet the requirements of both private and professional carers. They include a host of modern design innovations which ensure this lifter is lightweight, compact and efficient.

### Wangara 147JS Features & Specifications

- Rated capacity 147kgs
- Push handle height 950mm
- Low height (ground to hangar) 365mm
- Leg height 109mm
- Overnight battery charger
- Front castors 75mm twin
- Dual speed, hand held remote control
- Width (base internal) 470 - 960mm
- Tall height (ground to hangar) 1610mm
- Turning radius 1100mm
- Foot operated hydraulic leg spreader
- Rear castors 125mm single total locking

### Oxford Advance Folding Electric Hoist Features & Specifications

- Safe working load 227kg
- Max. overall length 139cm
- Minimum floor to 6pt spreader bar height 67.5cm
- Front twin castors 10cm
- No tools or special attachments needed to fold hoist
- Legs widest point is 108cm, narrowest is 67cm
- Advance's unique 'swan-neck' leg design allows the hoist to get even closer to the widest of obstacles
- Total assembly weight 31.7kg

**HIRE:** \$294\* for the first 4 weeks

**PURCHASE:** Oxford Advance \$1950



**GILLESPIES**  
HIRE & SALES SERVICE

"Comfort you deserve"

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027  
Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)