

LIFT & RECLINE ARMCHAIRS



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"



We are pleased to offer for HIRE and SALE a range of remote control, electrically operated 'Lift & Recline Armchairs' which offer the user a high level of comfort along with assistance to sit and stand.



Lift Function:

From the neutral, or normal sitting position, the simple remote hand control will raise the seat and tilt the chair forward to approximately 75% of a standing position. This is especially advantageous for those who have difficulty lifting their body weight from a seated to a standing position. The chair then remains in the 'lifted' position until the user returns to the chair, lowers themselves to the same position that they left the chair, then uses the remote control to lower themselves seated in the chair to the neutral seated position.



Recline Function:

The recline function operates in two stages. The first stage of movement, from the neutral seated position, moves the seat and backrest backward and down a little while extending the leg rest. This achieves a comfortable upright seated position with the legs elevated and higher than the hips. The second stage of movement lowers the backrest and raises the legrest simultaneously, to provide a reclined lying position many users find comfortable for rest or sleep.

HIRE: From \$294* for the first four weeks

PURCHASE: From \$1595

There are a variety of styles and fabrics available as well as variations in chair height and width. Six-way remote controls allowing individual operation of the legrest and backrest are also available as well as a range of options including castors for ease of movement; heat and massage units; neck comforters; battery power systems and fabric protection.

SPECIFICATIONS:

Please note that specifications do vary between models and actuating systems. Our most popular hire chairs have the following general specifications.

Seat Height	50cm
Seat Depth	47cm
Seat Width	50cm
Arm Height	17cm
Back Height	77cm
Lifting Capacity	140kg

*NB: Terms, conditions and deposits apply to hired goods

PHONE (02) 9411 2180

11-13 Elizabeth Street (PO Box 45) ARTARMON NSW 1570 AUSTRALIA



ambassador



Premier A1

Comfortable chair in traditional style built on heavy duty, single motor mechanism. This chair is upholstered in soft to touch fabric in four popular colours. Comfort and beauty blended together with strength and durability that will give you years of pleasure and peace of mind.

Features

- Rechargeable Battery Operating System
- Steel seat frame with no sag springs
- Two cushion backrest
- Integrated footrest overlay
- Pockets on both sides
- Headrest cushion
- Armrest covers

Options

- Hand control on left side at no cost.
- Integrated heat and massage system in lower lumbar (additional cost of \$400).
- Custom upholstery to suit your needs at additional cost based on individual order.

Warranty

12 months on lift & recline mechanism, all electrical components and structural timber frame.

Specifications

Lifting Capacity	up to 140 kg
Seat Height	50 cm
Seat Depth	48 cm
Seat Width	50 cm
Arm Height	17 cm
Back Height	77 cm
Overall Width	80 cm
Overall Height	115 cm
Overall Length	95 cm

Range of Colours



PA1ACO

PA1ASK

PA1ASA

PA1ALE

\$1,950

Note: Measurements are approximate only.

HIRE: \$294

FOR THE FIRST 4 WEEKS

Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

Assist-a-Lift Capri



POWERLIFT CHAIR

Made in USA

Assists in sitting and standing

3 Way mechanism

Easy to use hand controler

High quality Birch construction

Available in a range of colours

\$1,595

HIRE: \$294
FOR THE FIRST 4 WEEKS

Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE & SALES SERVICE

"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027
Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday

DAYS Stella Two Motor Riser Recliner Chair



days
Patterson Medical

Modern styling suitable for a wide range of homes
Dual motor enabling the user to raise their legs
without fully reclining the armchair
Available in Charcoal or Burgundy leather

\$2,695

Hire: \$334

FOR THE FIRST FOUR WEEKS

Dimensions:

Seat height: 47cm

Seat depth: 50cm

Seat width: 53cm

Overall width: 95cm

Overall depth (reclined): 172cm

Overall depth (upright): 90cm

Overall height: 120cm

Maximum user weight: 160kg



Queen

Ultra-comfortable chair with the convenience of reclining and standing
 Twin motor mechanisms operate backrest and footrest independently
 This chair is upholstered in soft to touch fabric in three popular colours

Features

- Battery Backup
- Two cushion backrest
- Integrated footrest overlay
- Pockets on both sides

Options

Hand control on left side at no cost.

Warranty

12 months on lift & recline mechanism, all electrical components and structural timber frame.

Specifications

Lifting Capacity	up to 110 kg
Seat Height	50 cm
Seat Depth	50 cm
Seat Width	50 cm
Arm Height	17 cm
Back Height	77 cm
Overall Width	80 cm
Overall Height	115 cm
Overall Length	95 cm

Note: Measurements are approximate only.



Range of Colours



QALE

QASA

QASK

\$1,950

HIRE: \$334

FOR THE FIRST 4 WEEKS

Terms, conditions and deposits apply to hired goods



TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au



Premier A4

Comfortable chair in modern style built on heavy duty, single motor mechanism. This chair is upholstered in soft to touch fabric in four popular colours. Comfort and beauty blended together with strength and durability that will give you years of pleasure and peace of mind.

Features

- Rechargeable Battery Operating System
- Steel seat frame with no sag springs
- Two cushion backrest
- Integrated footrest overlay
- Pockets on both sides
- Headrest cushion
- Armrest covers

Options

Hand control on left side at no cost.
Integrated heat and massage system in lower lumbar. Custom upholstery to suit your needs at additional cost based on individual order.

Warranty

12 months on lift & recline mechanism, all electrical components and structural timber frame.



Range of Colours



\$1,950

HIRE: \$294 For first 4 weeks

Terms, conditions and deposits apply to hired goods

Specifications

Lifting Capacity	up to 140 kg
Seat Height	45 cm
Seat Depth	48 cm
Seat Width	45 cm
Arm Height	17 cm
Back Height	77 cm
Overall Width	80 cm
Overall Height	115 cm
Overall Length	95 cm

Note: Measurements are approximate only.



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au



Premier A3

Comfortable bariatric chair in traditional style built on heavy duty, three motors mechanism. This chair is upholstered in soft to touch fabric in four popular colours. Comfort and beauty blended together with strength and durability that will give you years of pleasure and peace of mind.

Features

- Rechargeable Battery Operating System
- Steel seat frame with no sag springs
- Two cushion backrest
- Integrated footrest overlay
- Pockets on both sides
- Headrest cushion
- Armrest covers

Options

Hand control on left side at no cost.
Integrated heat and massage system in lower lumbar. Custom upholstery to suit your needs at additional cost based on individual order.

Warranty

12 months on lift & recline mechanism, all electrical components and structural timber frame.



Range of Colours ● ● ● ●



\$3,665

HIRE: POA

Terms, conditions and deposits apply to hired goods

Specifications

Lifting Capacity	up to 250 kg
Seat Height	50 cm
Seat Depth	50 cm
Seat Width	65 cm
Arm Height	17 cm
Back Height	77 cm
Overall Width	95 cm
Overall Height	115 cm
Overall Length	95 cm

Note: Measurements are approximate only.



TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday