



Stainless Steel Mobile Shower Commode



- Attendant-propelled transportation for showering or toileting
- Padded commode seat, footrest slides under frame to assist patient transfers
- Hand-made padded seat
- Durable stainless steel frame

QUALITY FEATURES

- Australian Made
- PVC slatted backrest & padded seat
- Swingback arms & sliding footrest
- Can be positioned for toileting
- Mounted on 4x125 mm castors, 2 front wheels are locking
- Pan carrier, stainless steel incl. (option)

SPECIFICATIONS

Armrest Height (from seat).....	240mm
Seat Height.....	590mm
Seat Width.....	445mm
Seat Depth.....	430mm
Back Height.....	365mm
Seat to Footplate	375mm
Overall Width	510mm
Weight Limit.....	120 kg

\$1,112

HIRE: \$204*

FOR THE FIRST 4 WEEKS

*Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE SALES SERVICE
"Comfort you deserve"



9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au