



# Mobile Shower Commode

Gillespies Hire & Sales Service  
is pleased to introduce the  
**'K-CARE Mobile Shower Commode**

With over 45 years industry experience we believe this model provides the best combination of features, high quality manufacturing and value for money available.



## QUALITY FEATURES

- Australian Made
- PVC slatted backrest & padded seat
- Swingback arms & sliding footrest
- Water resistant castors (2 x locking)
- Pre-zinc mild steel prevents rusting
- Can be positioned for toileting

Attendant propelled commode,  
shower chair and over toilet aid.

## SPECIFICATIONS

Armrest Height (from seat).....	240mm
Seat Height.....	590mm
Seat Width.....	445mm
Seat Depth.....	430mm
Back Height.....	365mm
Seat to Footplate .....	375mm
Overall Width .....	510mm
Weight Limit.....	120 kg

# \$734

**HIRE: \$140\***  
FOR THE FIRST 4 WEEKS

\*Terms, conditions and deposits apply to hired goods



**GILLESPIES**  
HIRE SALES SERVICE  
*"Comfort you deserve"*



# 9411 2180

11-13 Elizabeth Street, ARTARMON NSW 2064