

# Pride Lift Chairs

*Experience the Luxury*

## SPECIFICATIONS

|                             | C1                                    | C5                                    | C6   | C-101                     | LC-101 Euro Leather                                  |
|-----------------------------|---------------------------------------|---------------------------------------|--|---------------------------|--|
| <b>Back Width:</b>          | 51 cm (20")                           | 71 cm (28")                           | 71 cm (28")  | 72 cm (28.5")             | 58 cm (23")  |
| <b>Seat To Floor:</b>       | 41 - 46 cm (16" - 18")                | 46 - 51 cm (18 - 20")                 | 48 cm (19")  | 47 cm (18.5")             | 47 cm (18.5")  |
| <b>Seat Depth:</b>          | 51 cm (20")                           | 46 - 51 cm (18 - 20")                 | 53 cm (21")  | 51 cm (20")               | 51 cm (20")  |
| <b>Seat Width:</b>          | 46 cm (18")                           | 56 cm (22")                           | 56 cm (22")  | 51 cm (20")               | 51 cm (20") <sup>1</sup>                             |
| <b>Top of Back to Seat:</b> | 66 cm (26")                           | 66 - 71 cm (26 - 28")                 | 66 - 71 cm (26-28")  | 70 cm (27.5")             | 70 cm (27.5")  |
| <b>Weight Capacity:</b>     | 113 kg (250 lbs.)                     | 170 kg (375 lbs.)                     | 136 kg (300 lbs.)  | 147 kg (325 lbs.)         | 147 kg (325 lbs.)                                    |
| <b>Colour Options:</b>      | Arctic Blue,<br>Chocolate,<br>Oatmeal | Arctic Blue,<br>Chocolate,<br>Oatmeal | Arctic Blue,<br>Chocolate,<br>Oatmeal<br><b>Levante (vinyl):</b><br>Royal Blue | Cocoa,<br>Fern,<br>Sandal | <b>Leather:</b><br>Beige,<br>Brown,<br>Black,<br>Red |

  

|                             | D30   | 550                            | 560                           | 805               |
|-----------------------------|---|--------------------------------|-------------------------------|-------------------|
| <b>Back Width:</b>          | 74 cm (29")   | 89 cm (35")                    | 71 cm (28")                   | 76 cm (30")       |
| <b>Seat To Floor:</b>       | 48 cm (19")   | 48 cm (19")                    | 48 cm (19")                   | 48 cm (19")       |
| <b>Seat Depth:</b>          | 56 cm (22")   | 48 cm (19")                    | 48 cm (19")                   | 51 cm (20")       |
| <b>Seat Width:</b>          | 56 cm (22")   | 53 cm (21")                    | 48 cm (19")                   | 51 cm (20")       |
| <b>Top of Back to Seat:</b> | 66 cm (26")   | 68.5 cm (27")                  | 63.5 cm (25")                 | 69 cm (27")       |
| <b>Weight Capacity:</b>     | 170 kg (375 lbs.)   | 170 kg (375 lbs.)              | 170 kg (375 lbs.)             | 136 kg (300 lbs.) |
| <b>Colour Options:</b>      | Ocean,<br>Burgundy,<br>Spruce,<br>Latte<br><b>Levante (vinyl):</b><br>Burgundy,<br>Royal Blue,<br>Pine Mountain | Almond,<br>Merlot,<br>Rivergum | Ocean,<br>Burgundy,<br>Spruce | Blue,<br>Expresso |

## BENEFITS

- Engineered furniture laminate/hardwood frame (C1, C5 and C6 all have steel frames)
- Patented, quiet & smooth lift system
- Large side pocket
- Hand control can be easily moved to left hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup
- 2-year limited warranty

Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

21 Healey Rd.  
Dandenong, Victoria 3175

Tel: 03 9706 4611  
Fax: 03 9706 4622  
www.pridemobility.com



**Pride**  
Mobility Products  
Australia Pty. Ltd.  
— LIVE YOUR BEST! —

**PRIDE**  
— LIVE YOUR BEST! —



# Pride® Lift Chairs

*Experience the Luxury*

# Pride Lift Chairs

*Experience the Luxury*

Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any décor, Pride Lift Chairs also feature a 2-year limited warranty. The state-of-the-art lift system utilized in Pride Lift Chairs will recline you further back than any other lift chair. Plus, Pride Lift Chairs are equipped with an emergency battery backup system should the power fail in your home.

## LIFT CHAIR DESCRIPTIONS

### C1

3-position, petite to small lift chair with fully removable upholstery components

### D30

3-position, medium to large chaise lounge

### C5

3-position, medium to large chaise lounge with fully removable upholstery components

### 550

3-position, medium to large chaise lounge with heat and massage

### C6

Infinite position, medium to large chaise lounge with fully removable upholstery components and dual motors which enable independent recline/footrest function

### 560

3-position, small to medium chaise lounge

### C-101

3-position, medium to large lift chair

### 805

2-position, medium to large wall hugger chaise lounge which requires only 10 cm (4") of space between wall and chair

### LC-101 Euro Leather

3-position, medium to large lift chair

## AVAILABLE FABRICS <sup>1</sup>



Spruce Ocean Rivergum



Burgundy Latte Almond



Blue Merlot Oatmeal

### SUEDE



Chocolate Arctic Blue



Fern Cocoa



Sandal Espresso (Padded Suede)

### LEATHER



Black Beige



Red Brown

### LEVANTE VINYL



Burgundy Royal Blue



Pine Mountain

## LIFT CHAIR MODELS



### C1

Split-back design  
Shown in Arctic Blue



### C5

Sewn Pillow-back Design  
Shown in Arctic Blue



### C6

Sewn Pillow-back Design  
Shown in Chocolate



### C-101

Button-back Design  
Shown in Sandal



### LC-101 Euro Leather

Split-back Design  
Shown in Beige Leather



### D30

Split-back Design  
Shown in Burgundy



### 550

Split-back design  
Shown in Almond



### 560

Pillow-back Design  
Shown in Spruce



### 805

Sewn Split-back Design  
Shown in Blue

**PLEASE NOTE:** 2-position is a TV recline.  
3-position is a full recline.  
Infinite-position is a full recline with independent legrest function.



The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

<sup>1</sup> Not all fabrics are available for every model. Please see Specifications for details.