

SEATING



Day chair - High Back with Air Cushion

Features:

High back chair with air cushion seating
Adjustable legs
Lumbar support
Weight capacity 152kgs

Specifications:

Seat width 52cm
Seat depth 46cm
Backrest height 70cm
Adjustable seat height 41-56cm

HIRE: \$154*
for the first 4 weeks

PURCHASE: \$569



Carequip 'Hunter' Chair

Features:

Low back orthopaedic chair
Ideal for dining use or household activities
Adjustable seat height and angle provide optimum hip support
Curved padded backrest provides postural support and comfort
Comfortable padded armrests assists in standing and transfers
Powder coated steel frame
Maximum user weight: 100kg

Specifications:

Seat: width: 440mm
Seat depth: 440mm
Seat height: 495-645mm
Arm height: 220mm
Back height: 450mm
Overall width: 580mm

HIRE: \$94*
for the first 4 weeks

PURCHASE: \$299

The 'Revolution' Chair

The revolution chair incorporates a unique, patented "swivel-slide" technology allowing easy access to the chair and to a table edge without moving its legs.

The lock/release mechanism allows the seat to gently slide back from the table so it can swivel left or right for easy exit.

Specifications:

Seat width: 450mm Overall width: 505mm
Seat height: 480mm Overall depth: 560mm
Seat depth: 485mm
Maximum user weight: 120kg

PURCHASE: \$725



Made in Australia



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au