

# Sterling Sapphire<sup>2</sup> Scooter



## Features:

- Full size scooter than comes apart in seconds like a smaller boot scooter to offer you the best of both worlds
- Infinitely adjustable tiller
- Wrap-around handlebars
- Finger-thumb directional control
- Adjustable speed up to 6.4kph
- 150kg safe working load

## Purchase

# \$1,950

## HIRE: \$490\*

FOR THE FIRST 4 WEEKS

\*Terms, conditions and deposits apply to hired goods

## Specifications:

- Length: 121cms
- Rear width: 59cms
- Battery Capacity: 35Ah
- Max Safe Slope: 14 degrees at 113kg  
10 degrees at 150kg
- Turn Radius: 70cms
- Seat - Base Height: 50-62cms
- Wheel Diameter Front: 22.5cms
- Wheel Diameter Rear: 26.5cms
- Charger Off-Board: 5Amps
- Maximum speed: 6.4kph
- Overall Weight: 79kg
- Controller (PGDT): S-Drive
- Ground Clearance: 9.5cms



**GILLESPIES**  
HIRE & SALES SERVICE  
"Comfort you deserve"

**TEL: 02 9411 2180**

Email: [sales@ghss.com.au](mailto:sales@ghss.com.au)

Fax: (02) 9411 2173

Web: [www.ghss.com.au](http://www.ghss.com.au)