



**GILLESPIES**  
**HIRE SALES SERVICE**

*"Comfort you deserve"*

## Stocking Aid

A stocking aid which will apply support hosiery quickly and evenly. Grooves catch the compression stocking and prevent it from popping off. Works like a traditional sock aid - just pull the stocking over the shell, insert the foot and pull the cord to bring the hose up onto the leg. Suitable for Class 1 and 2 stockings. Recommended for below knee stockings but may still be suitable for some above knee models.

Small: 130mm wide x 140mm high (5 1/4 x 5 1/2").

Large: 150mm wide x 145mm high (6 x 5 3/4").

This product may be referenced by 09-105-9971, 091059971 or A75441  
**STOCKING AID ROLYAN COMPRESSION LARGE**



**\$105.95**

Visit our website for more details - [www.ghss.com.au](http://www.ghss.com.au)

**PHONE (02) 9411 2180**

Open Monday to Friday 9am - 5.30pm & Saturday 9am - 12 Noon  
13 Elizabeth Street, Artarmon NSW 2064 Australia