

PORTABLE GRAB BAR

TELESCOPING GRIP



MARKET FACTS

Falls are the most common accident around the home

More than a third of adults aged 65 years or older fall each year*

Of those who fall, 20% to 30% suffer moderate to severe injuries that reduce mobility and independence, and increase the risk of premature death*

A recent study of people aged 72 and older found that the average health care cost of a fall injury was \$19,440*

The average installation and materials cost for installing a permanent grab bar with blocking is \$537**

* CDC **Industry Labor Council

PRODUCT SPECIFICATIONS

Item Description: Telescoping Grip
Portable Grab Bar

Color: White

Length: Medium Length 19" -24"
Part Number: BC500M
UPC: 898387001013

Large Length 24" - 29"
Part Number: BC500L
UPC: 898387001020

Suction Pad Diameter: 4.7"

Shipping Weight: 3.75 lbs

Patent Pending



FEATURES

Installs and is removable without tools or professional installation

Durable Commercial Grade Suction Cups

Industry leading superior functionality

Contemporary stylish design

Will not damage or deface property

Adjustable length offers adaptability to most any tub/shower tile configuration

Telescoping handle allows for a much longer hand holding area

Offered in two sizes with adjustable overall lengths ranging from 19"-29".

Limited lifetime warranty

Because it is portable and will not damage or deface property, it is great for:

- Bathtub and Showers
- Vacations
- Hotel Stays
- Winter Destinations (perfect for snow birds)
- Renters
- Rehabilitation
- Anytime a permanent grab bar is not possible.



TEL: 02 9411 2180 WEB: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday