

TOILET AIDS

Over Toilet Aid



Tubular Steel fabrication
Zinc treated frame to resist rusting
Adjustable Seat Height 480 - 610mm
Width - Arms 450mm / Base 540mm
Weight capacity 100kgs
Armrests to assist sitting & rising
Also in 'Wide/Heavy Duty' version

HIRE:
from \$55* for the first 4 weeks

PURCHASE:
K-Care **\$135**

Toilet Surround



Tubular Steel fabrication
Phosphate treated to resist rusting
Adjustable Height 660 - 780mm
Width between armrests 450mm
Armrests to assist transfers

HIRE:
from \$50* for the first 4 weeks

PURCHASE:
K-Care **\$110**

Raised Toilet Seat



Moulded plastic seat and lid
Available in 50mm, 100mm & 150mm
Clamps to both sides of toilet pan
Strong & lightweight
Weight capacity 180kgs
Simple to install and remove
Fits most toilet suites

HIRE: from \$50*
for the first 4 weeks

PURCHASE:
100mm **\$115.00**

Commodes

For persons who need a toilet facility placed more conveniently, to avoid transfer to the bathroom, we hire and sell a range of 'Commodes' which provide a toilet facility able to be positioned next to the bed or in any level area. Our Commodes are supplied with a 'pot & lid' or 'bed pan' dependant upon the style.



Please refer to, or request, our standard 'Commode Brochure' for a more detailed description of both the equipment and associated charges.

GRAB RAILS:

We stock a range of steel grab rails which may be fitted to stable surfaces in the bathroom (walls & cabinets etc) to give support and aid stability. We supply specialised stainless steel rails, 'made to order' in various shapes to suit your particular requirements.

ACCESSORIES:

Also available are a wide variety of toilet & incontinence products including;

- Urinals (Male & Female)
- Bed Pans & Slipper Pans
- Mattress protectors (S,D,Q)
- Absorbant Bed Pads (Kylie / Drycare)
- Absorbant Chair Pads
- Blueys & Maxi Blueys

Please feel free to call and discuss your specific needs with our friendly staff.
We look forward to being of service.



*NB: Terms, conditions and deposits apply to hired goods

PHONE (02) 9411 2180

11-13 Elizabeth Street (PO Box 45) ARTARMON NSW 1570 AUSTRALIA