



'Vito Plus' Transit Wheelchair

202



- Swingaway removable legrests and fixed armrests
- Folding backrest and frame make a compact unit for storage & transportation
- Attendant hand brakes provide more control for the carer
- Desk armrests allow wheelchair to be pushed close to desk or table
- Adjustable nylon buckle seat belt for additional patient safety
- Heel loops
- Breathable nylon upholstery
- Silver paint finish
- Quality lightweight aluminium

Seat Width:	460mm
Seat Depth:	400mm
Seat Height:	480mm
Back Height:	460mm
Arm Height:	210mm
Push Handle Height:	910mm
Overall Width:	530mm
Gross Weight:	11.7kg
Net Weight:	10.45kg
Max User Weight:	110kg