



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027
Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday



DV8 Steerable Knee Walker

790 Steerable Knee Walker, 1/cs

FEATURES:

- Knee walker can be steered for increased maneuverability.
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot.
- Deluxe lever brakes ensure safety.
- Tool free height adjustment.
- 8" casters are ideal for indoor/outdoor use.
- Limited Lifetime Warranty.
- Lifetime Warranty on brake cable.

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Depth	Pad Height (highest)	Pad Height (lowest)	Pad Width	Pad Depth	Product Weight	Shipping Weight	Weight Limit
790	37.79"	33.07"	16.9"	31.49"	21.81"	18.81"	7.08"	14.17"	22.70 lbs.	28.88 lbs.	300 lbs.