

HOSPITAL BEDS & PATIENT LIFTERS

We hire and sell a range of bedroom equipment for both temporary & permanent requirements, such as remote controlled, electrically operated hospital beds and patient lifters to suit most needs. We also carry a range of pressure care overlays (please refer to our 'Pressure Care' brochure), and alternating air mattresses for those requiring additional relief. Please feel welcome to call and discuss your specific needs with our experienced staff.

HOSPITAL BEDS

K-Dee Standard Bed



K-Dee Standard Bed Features & Specifications

- 4 section bed
- Easy to assemble
- Can be folded for easy transport
- Adjustable back and foot rests (back 70°, foot 26°)
- Tilting
- Knee break
- Four locking castor wheels (4"/10cm)
- Linak Motors (Linak actuators)
- Dual height adjustment (33-73cm or 40-80cm)
- Maximum user weight 180kg

HIRE: \$540*
for the first 4 weeks

PURCHASE: \$2,220 Code: B014R
\$650 mattress

PATIENT LIFTERS

We hire electric patient lifters with slings to aid carers in the safe and convenient transfer of mobility impaired persons. The patient lifters are designed to meet the requirements of both private and professional carers. They include a host of modern design innovations which ensure this lifter is lightweight, compact and efficient.

Oxford Advance Folding Electric Hoist Features & Specifications



- Safe working load 227kg
- Max. overall length 139cm
- Minimum floor to 6pt spreader bar height 67.5cm
- Front twin castors 10cm
- No tools or special attachments needed to fold hoist
- Legs widest point is 108cm, narrowest is 67cm
- Advance's unique 'swan-neck' leg design allows the hoist to get even closer to the widest of obstacles
- Total assembly weight 31.7kg

HIRE: \$294* for the first 4 weeks

PURCHASE: Oxford Advance \$2,590 Code: OXSADV



GILLESPIES
HIRE & SALES SERVICE

"Comfort you deserve"

TEL: 02 9411 2180


Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027
Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday

Classic FS Bed



'DAYS' Classic FS Bed

FEATURES

- 4-section mattress support, 2-motor adjustment, easy transport
- Client position assistance pole with triangle
- Wooden head end, foot boards, side rails in beech finish
- Handset includes locking facility to maintain position
- Removable mattress retainers at head and foot sections
- Ergonomic adjustment of head and foot sections
- Additional positioning of foot section by ratchet fittings
- Emergency lowering by quick-action bolt mechanism
- Transport lock
- Dual action headrest and knee break adjustment on handset, to help reduce shear when profiling

TECHNICAL DATA

- Total weight....(4-piece) 91.5kg
- Heaviest part....20.1kg
- Outer dimensions W/L....102.3 x 212.5cm
- Max safe working load....175kg
- Max patient weight....140kg
- Height adjustment....40cm to 80cm
- Adjustment angle of back section....0° - 70°
- Adjustment angle of foot section....0° - 20°

QUALITY

- EN 1970:2000 | DIN EN 1970:2000-12
- EN 60601-2-38:1996 + A1:2000 (as basis)
- DIN EN 60601-2-38:2001 (as basis)
- CE 93/42/EEC on medical devices

Code: DAYSFSBEDSULTRA

*Great value,
Fast delivery!*



HIRE: \$540*

FOR THE FIRST 4 WEEKS
Hire cost includes mattress,
side rails and self-help pole

\$2,200

SALE PRICE MATTRESS NOT INCLUDED

*Terms, conditions and deposits apply to hired goods



GILLESPIES
HIRE & SALES SERVICE
"Comfort you deserve"

TEL: 02 9411 2180

Email: sales@ghss.com.au

Fax: (02) 9411 2173

Web: www.ghss.com.au

11-13 Elizabeth Street, (PO Box 45) Artarmon NSW 1570 Australia - GHSS Pty Ltd ABN 94 082 770 027

Open Monday to Friday 9am - 5.30pm & Saturdays 9am - 12 Midday